Türkçe / English

KK-01 09.2022/Rev:06

KULLANIM KILAVUZU

Instructions for Use

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Store in a dry place

Protect from sunlight

CE marking according to Directive 93/42 EEC

Manufacturer

Production Date: Located on the package.

SN Serial Number: Located on both the product and the package.

REF Catalog Number: Located on both the product and the package.

INTENDED USE AND INDICATIONS

Reusable Surgical instruments means an intrumen intended for surgical use in cutting, drilling, sawing, scratching, scraping, clamping, retracting, clpping or similar procedures without a connection to an active device and which is intended by the manufacter to be reused after appropriate procedures such as cleaning, disinfection and sterilisation have been carried out.

CONTRAINDICATIONS

There no known contraindications however contraindicated for all other applications except their intended use! The conditions that can risk the success of the operation such as injuries or accompanying diseases should be evaluated by the user. Risks should be taken into account for anatomical structures in the planned area.

GENERAL SAFETY WARNINGS

- Read, the instructions for use before each use.
- Keep this instructions for use throughout the product use.

- Products are **NON-STERILE.** It must be cleaned and sterilized before use according to the instructions below.

- Used only by people the requisite training, knowledge and experience.
- Do not use out of its intended use.
 Before each use, instruments should be checked / inspected. See, Inspection & Check
- Do not use out of sight or in situations where vision is not available.
- In cases with infectious disease risk, the legislation determined by the relevant authority institution must be followed.
- Before each use, check the intercompatibility of instruments that are used with another device.
- Keep the HF instruments away from explosives and flammables.

WARNINGS REGARDING CLEANING, DISINFECTION AND STERILIZATION

- For the cleaning, disinfection and sterilization of products, national and international standards, directives and clinical hygiene instructions must be followed.
- Use cleaning and disinfectant agents according to the manufacturer's instructions.
- Apply cleaning processes (manual or mechanical) to newly purchased instruments before the first sterilization process.
- For patients with a diagnosis or suspicion of Creutzfeldt-Jakob disease (CJD) or variants, national regulations regarding disinfection and sterilization of products must be followed.
- Handle with validated methods and calibrated autoclaves and machines at all times.
- Use Personal Protective Equipment (PPE) when working with contaminated instruments.PPE: apron, mask, glasses, gloves, shoe covers....
- Prefer mechanical cleaning cleaning instead of manual cleaning as it gives more reliable results.
- Handle with care. Do not drop or apply excessive forces to the instruments during reprocessing.

Cleaning and Maintenance Products / Tools

warning **The use of metal brushes or clothes cause corrosion !**

warning **The use of chlorine solutions (Saline solutions, bleach, Ringer's solution) cause pitting corrosion or stress corrosion cracking !**

- Use cleaning and disinfection agents appropriate for aluminum, high-grade steel.
- Use of oxidizing chemicals (eg. H2O2) with PVD coated instruments causes fading and chemical stripping.
- Do not use hydrogen peroxide solutions (H2O2) on instumens with titanium and its alloys..
- Long-term processes with disinfectants based on alklamine active ingredient in silicone elastomers (silicone-coated products) may cause hardening
- Use cleaning agents with approved efficacy (VAH or FDA approval or CE mark)
- Cleaners with a disinfectant effect (manual cleaning) must be in accordance with EN 14885 or equivalent national guidelines.
- Use soft, lint-free fabrics, paper cloth, plastic, soft bristle brushes and water spray guns
- Do not use aldehyde based cleaning products due to its protein fixing properties.
- Use non-flaming cleaning agents and solutions for HF instruments
- Cleaning solutions should be changed on daily basis. If the solution gets contamined visibly, it can be changed earlier.
- If powder products are used, the powder should be completely dissolved in water before use. Undissolved particles may cause changes in instrument surfaces and clog narrow tube instruments.
- Compressed air drying is more effective than drying with fabric or cloth.

Ultrasonic cleaning is recommended to contribute to both manual and mechanical cleaning especially for lumened and articulated instruments

- Suitable detergent should be selected for ultrasonic cleaning. (Non-foaming)
- Separate stainless-steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis.
- Be sure that products are free of cleaning agent at the final rinse of ultrasonic cleaning step.
- Micro surgical instruments must be places on suitable racks in order to prevent damage.

Water qualifications

Softened water*: Water hardening substances (calcium and magnesium cations) are reduced.

Demineralised water*: Mineral substances largely removed by reverse osmosis or cation & anion exchangers or electrode ionization or distillation.

For steam sterilization and the last rinsing of cleaning processing, limit values for feed water quality is specified in EN 285 and ISO 17665.

Non-conforming water properties cause:

- Scaling, lime deposits,
- Brown-red deposits, secondary rust,
- White-grey, colored appearance, thin scaling,
- Pitting,
- Spotting and scaling.

*See www.ak-i.org, Publications, Red Brochure for water properties (softened water, demineralized water)

CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

Point of Use

- Remove excess biologic soil from the instruments with a disposable wipe.
- Rinse non-visible surfaces with demineralised water, with a disposable syringe.
- Place dry products in sealed waste containers.

Containment and Transportation

warning Wet disposal conditions increase the risk corrosion and microbial growth. The residues of cleaning agents interfere with the following processes.

- Transport the instruments to decontamination are in closed containers at dry state.
- Keeping the products at a dirty state more than 6 hours may cause corrosion and fail cleaning
 processes. Used instruments should be processed within 6 hours at dry disposal.
- Prevent slipping of the instruments by using special carrying tools during transportation.

Preparation for Cleaning / Precleaning

warning

Do not exceed 25 °C for pre-cleaning temperature.

- Rinse the instruments to remove coarse contaminations.
- Disassemble the detachable instruments to preparing for cleaning.
- Ultrasonic washing is recommended before disinfection to increase the cleaning efficiency of heavily soiled or difficult to clean instruments such as lumened/tubular, micro and articulated instruments for mechanical support.
- Pre-treating of HF instruments with hydrogen peroxyde solution is recommended.
- If a pre-cleaning step is necessary before automatic cleaning / disinfection, follow the 1st and 2nd process steps in the "Manual Cleaning / Disinfection" section of this document. Then send the instruments to the automatic washing-disinfection process by dry transfer.

Manual Cleaning / Disinfection

warning

Do not exceed the maximum permitted cleaning temperature of 45 °C. Otherwise protein denaturation may ocur.

warning

Use a non-protein fixing cleaning solution does not contain alcohol, aldehyde, phenol and QUAT with a disinfection effect according to the manufacturer's instructions.

warning

All scrubbing should be performed below the surface of the solution. Do not apply manuel cleaning to the cutting instruments or sharp pointed instruments to avoid physical demage to the operator !

	Process Step	Water Quality	Heat	Time
1	Cleaning/disinfecting[1]	Softened water*	Room Temperature	15 min
2	First Rinsing	Softened water*	Room Temperature	1 min
3	Disinfecting[1]	Softened water*	Room Temperature	15 min
4	Final Rinsing	Demineralized Water*	Room Temperature	1 min
5	Drying	-	Room Temperature	**

1. Pour a 25 ml dose cleaning/disinfectant solution in 5 liters of water (%0,5 dilution). Completely immerse instruments in the solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Flush the solution into non-visible surfaces and holes of lumened and tubular instruments such as cannulas with a disposable syringe (20 ml). While soaking scrub surfaces using a suitable cleaning brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to hinged joints, box locks, instrument teeth, rough surfaces and areas with moving components or springs. Lumens, blind holes and cannulas should be cleaned using a snug fitting round suitable brush. Insert the brush into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.

 Note: If the instrument is too large (such as table-mounted bars) to be immersed into the solution, a proper disinfectant cleaner should be sprayed out surface of the instrument.

2. Remove the instruments from the solution and rinse with softened water for a minimum of 1 minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas. Allow excess water to drip of from the instrument prior to disinfectant solution.

3. Pour a 125 ml disinfectant solution dose in 5 liters of water (%2,5 dilution). Immerse fully in the disinfectant solution. Move the joints at least 5 times beneath the solution. Flush the solution into lumens with a disposable syringe (20 ml).

4. Rinse with running demineralized water for at least 10 seconds close and open the joints at least 5 times, flush lumens with a disposable syringe. We recommend using fully desalinated water (demineralized, drinking water quality at least) to avoid water spots.

5. Dry instruments with a clean, absorbent non-shedding lint free cloth or compressed air to remove moisture from lumens, holes, cannulas and difficult to access areas.

[1] Anios - Aniosyme DD1 Solution was used to verify cleaning.

- Check for cleanliness visually at the end the cleaning process especially the hinges, lumens, ratchets, teeth and etc. Repeat the cleaning process if necessary.
- Apply maintenance products (see **Maintenance & Care**)
- Control visually and functioning of the instruments according to the **"Inspection & Check"** part of this document.

Mechanical cleaning and thermal disinfecting

warning

Use low alkaline cleaners (8.0-11.0 pH). Do **not** use high alkaline value (>11 pH) cleaning agents.

warning

The washer disinfecfer must have either FDA approval or CE mark (validated according to EN ISO 15883)

- Place instruments in a washer/disinfector. Follow the washer / disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure; e.g. open all instruments, place concave instruments on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfector is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer's instructions. Microsurgical instruments must be stored on special holders or suitable racks in order to prevent damage.
- Arrange the position of the instruments in the washer/disinfector to prevent mechanical damage through contact.
- **Avoid overloading** to make sure the surface of the instrument effectively gets in contact with cleaning solution.
- To prevent corrosion, unload the washer/disinfector immediately when the cleaning program ended.
- The following minimum wash cycle parameters are validated and recommended:

	Process Step	Water Quality	Heat	Time
1	Pre-wash	Softened water*	<25 °C	3 min
2	Cleaning	Demineralized Water*	40 – 55 °C	10 min
3	First rinse / Neutralization	Demineralized Water*	>10 °C	2 min
4	Intermediate rinse	Demineralized Water*	>10 °C	1 min
5	Thermal disinfection / Final rinse	Demineralized Water*	90 °C	5 min
6	Drying	-	-	**

- Check for cleanliness visually at the end the cleaning process especially the hinges, lumens, ratchets, teeth and etc. Repeat the cleaning process if necessary.
- Apply maintenance products (see **Maintenance & Care**)
- Control visually and functioning of the instruments according to the **"Inspection & Check"** part of this document.

* See "Water qualifications" section of this document for water properties.
** Dry time is dependent upon the load size placed into the washer/disinfector.

Maintenance & Care

warning

Be sure to apply care products (lubricant) to metal sliding surfaces, moving parts, joints/box locks prior to **"Inspection & Check"**

warning

Raw meterials used in care products must comply with pharmacopeias such as USP/NF, Ph. Eur, and national and international cosmetics regulations. **Do not use** silicone oil based products!

- Ensure that the instruments are cleaned effectively before maintenance. Visually and physically inspect the instruments for blood and protein soils. It is recommended to verify the cleaning with a chemical protein test, especially on surfaces that are difficult to clean or invisible such as ratchets, hinges, teeth, lumens.

 Care products must be:

 - Paraffin/white oil based,
 - Biocompatible,
 - Suitable for steam sterilization and vapor-permeable.

- Allow the products to cool down to room temperature.
- Assemble the diassembled parts. During the assembly for detachable rongeurs/punches make sure marked numbers are matched.
- Apply the care product directly to joints, threads and friction surfaces by using a squirt oiler, an oil pen, or through targeted application of drops of oil
- Actuate moving mechanisms a few times to distribute the care product evenly.
- Wipe the excessive care product with a lint-free cloth.

warning

Failure to apply sufficient care products to moving parts causes friction related damage and corrosion!

Inspection & Check

Verify the following:

- Check the laser marking of the product. It must be easily readable visually or mechanically.
- There must not be any crack, fracture, wearing, burr, discoloration, spotting, rust or corrosion on the surface of the product.
- There must not be any loose, bent, broken, cracked, worn, or fractured components.
- There must not be any discoloration or deterioration in coating integrity on coated, painted and isolated coated products.
- Lumened/tubular instruments must be checked against blockages with a suitable cleaning brush. If a blockage is detected, reprocessing must be repeated.
- Cutting and puncturing instruments should not be dull.
- Moving parts must be able to move freely without getting stuck or rubbing.
- Products with spring-loaded components should regain their original form after use.
- Matching parts must be compatible with each other for the assembled products.

warning

If any of the above controls are not confirmed, do not use the product, contact with Bahadır Medical Devices (manufacturer).
Damaged products must **only** be maintained by the manufacturer.

Packaging prior to Sterilization

- Pack the products according to the determined sterile packaging. *
- Lock the products with ratchet locks in the first notch.
- Ensure that all cutting edges and fine tips are protected.

***The Bahadır Sterilization Containers were used in the validated sterilization method.**

Sterilization Method

Moist heat/steam sterilization is the preferred and recommended method for Bahadır Reusable Surgical Instruments.

- The steam sterilizer must comply with the EN 285 standard and be verified according to the EN ISO 17665 standard.
- Ensure that the maximum permitted load according to the manufacturer's instructions is not exceeded.
- For steam sterilization, limit values for feed water quality as specified in EN 285 and ISO 17665.
- Ethylene oxide or other sterilization methods are **not** recommended.
- The use of chemical indicators causes stains on the instrument surface.

Sterilization Method	Heat	Time
Prevacuum	134 °C[1,3]	5 min[2,3]

[1] Validated exposure time required to achieve a 10-6 sterility assurance level (SAL)
[2] Validated exposure temperature required to achieve a 10-6 sterility assurance level (SAL).
[3] Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in this table.

Storage

- Storage and distribution of the product should be maintained by the designated staff according to the rules of the instutition which is used.
- Newly purchased products should be stored in a dust and moist-free environment.
- Surgical instruments should **not be** stored near to chemicals.
- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.

Life Span

- As long as the instruments are regularly maintained, used according to their intended use, storage and maintenance conditions are followed, they can serve for 7 years depending on the average usage intensity.
- Products which can not pass from functional and visual inspection can not be used.

Disposal

- The disposal of a product that comes into contact with the human body and its packaging should be done in accordance with the national laws of the country where the product is used.

Technical Service

- Guarantee term: 2 years
- The user is responsible for any adverse situation that may arise from the use of products whose warranty expires.

warning

In case of any adverse event related to the device, Bahadır Tıbbi Aletler AŞ (manufacturer) or local authorities should be reported.

Technical Service Contact Details

Manufacturer
BAHADIR TIBBİ ALET CİHAZ VE İNŞAAT MAKİNA SANAYİ VE TİCARET ANONİM ŞİRKETİ
Şabanoğlu OSB Mah. Vali M.Erdoğan Cebeci Bul. No:29/1
Tekkeköy / Samsun

Tel: +90 362 366 70 00
Fax: +90 362 366 70 05
info@bahadir.com.tr / www.bahadir.com.tr

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