Türkçe / English

KK-02 05.2023 / Rev:05

KULLANIM KILAVUZU

Instructions for Use

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ENG Sterilization Container Systems

Symbol Explanations

Read instructions before use!		CE marking according to Regulation (EU) 2017/745	
Warning!		Manufacturer	
Non-Sterile		Production Date: Located on the package.	
Do not reuse		Serial Number: Located on both the product and the package.	SN
Store in a dry place		Lot/Batch Number: Located on both the product and the package.	LOT
Protect from sunlight		Catalog Number: Located on both the product and the package.	REF

1. GENERAL SAFETY WARNINGS

- Read, the instructions for use before each use.
- Keep this instructions for use throughout the product use.

- Products are **NON-STERILE!** It must be cleaned and sterilized before use according to the instructions below.

- The products must be used only by people with the requisite training, knowledge and experience.
- Do not use other than intended use!
- Before each use, products should be checked/inspected. (See **Inspection & Check**)
- In cases with infectious disease risk, the legislation determined by the relevant authority.

2. INTENDED USE AND INDICATIONS

The Bahadır Sterilization Containers are used for the reprocessing of medical devices, packaging, transportation, and storage of instruments prior to, during, and after sterilization.

Indicated for 3x fractionated pre-vacuum procedures.

3. CONTRAINDICATIONS

- All other applications are contraindicated except for the steam sterilization methods.

Comparison table to demonstrate the suitability of container types according to sterilization methods:


Container Types Lid / Bottom	Sterilization Methods				
	Pre-Vacuum	Gravity	Fractional Vacuum	Fractional Flow	Steam Injection
Perforated / Perforated	✓	✓	✓	✓	✓
Perforated / Nonperforated	✓	-	✓	✓	✓
Nonperforated / Perforated	✓	-	✓	✓	✓
BioBarrier Systems	✓	-	✓	✓	✓
Nonperforated / Nonperforated *	-	-	-	-	-

* Removal containers (without filters and holes) are contraindicated for sterilization of medical devices. It can only be used for the transport and disposal of medical waste. If it requires to be sterilized, place it in the sterilizer in the open position. If not, the forms may be deformed due to vacuum!

4. TECHNICAL DESCRIPTIONS

- Bahadır Container Systems are manufactured from anodized aluminum alloy and surgical stainless steel according to EN 868-8, DIN 58952-2, and DIN 58952-3 standards.
- Bahadır Container Systems has been tested with tests in accordance with EN ISO 868-8 Annex-D standard and is produced in a design that can be stacked on the shelf and in the sterilizer.
- Bahadır Sterilization Container Systems consist of containers, filters, wire/sheet metal trays, and accessories. Accessories are a component of the container and are not used externally.
- Bahadır Container Systems has different model and size variants. Perforated/non-perforated versions and different filter options are available depending on the sterilization method and preferences. Container lids can be preferred with safety lids that protect the filter against damage and reduce the risk of contamination. In addition to versions that can be disassembled by the user, systems integrated into the lid are also available.

4.1 Structure


(1) Safety Lid*
(2) Lid
(3) Filter/Filter Retainer**
(4) Container (Bottom)
(5) Front Panel
(6) Indicator Slot
(7) Label Slot
(8) Handle

* (see 5.2 Safety Lids)
** (see 5.1 Filters)

5. USAGE / APPLICATION


- Use The Bahadır Containers, together with the accessories and other components in this catalog (CT20). Please review this catalog for the suitability of the trays to be used in the containers.

5.1 Filters


- **(1) Paper Filter:** The filter is disposable. It is attached on the lid with the filter retainer.
- **(2) PTFE Filter:** The filter has 1500 cycles of use. It is attached on the lid with the filter retainer.
- **(3) BioBarrier Filter:** The filter is reusable. It is integrated on the lid.
 - The BioBarrier filter upper cover is removed by turning counterclockwise, applying pressure from the center point. It is fixed by turning clockwise by applying pressure from the center point for locking.
- **(4) Filter Retainer:** For the usage of paper and PTFE filter.
 - Release the filter retainer lock with the help of the button shown on the image. After the filter is placed in the container, place the filter retainer on it. Perform the locking operation by applying pressure from the center point. After the hear click sound, be sure that the locking by checking manually.

5.2 Safety Lids

- **(1) Demountable Safety Lid:** Suitable for container lids with disposable filters.
 - Release the lock mechanism by pressing the button as shown in the image. Lift the lid to remove it. When refitting, start with the non-locking edge, then push the other edge onto the lid.
- **(1) Integrated Safety Lid:** It comes as standard in W Series lids with reusable filters and M Series lids, securely fastened to the lid with a central screw.

5.3 Front Panel


- **(1) Indicator:** When the temperature reaches 134 °C during sterilization, it is verified that the appropriate temperature has been reached by the dark brown color change in the chemical indicators. Check the color change on the indicator after each sterilization process.
 - There is also a chemical indicator on the disposable paper filter, perform internal control with that during the aseptic presentation.
- **(2) Security Seal:** The security seal ensures the safety of the sterile package against unauthorized opening. It is placed to the locking hinge before sterilization as shown in the image and opened by break off.
- **(3) Labels:** The labels are used to categorize information of the surgical set, department, etc.

Do not reuse

Do not reuse disposable components such as paper filters, indicators and security seals!

5.4 Silicone Gaskets

- Silicone gaskets are seated on the retaining groove and filter retainer and ensure microbiological sealing.
 - The gaskets should only be replaced by the manufacturer and authorized personnel.

⚠ 6. WARNINGS REGARDING CLEANING, DISINFECTION AND STERILIZATION

- For the cleaning, disinfection and sterilization of products, national and international standards, directives and clinical hygiene instructions must be followed.
- The containers must be properly cleaned and sterilized prior to first use.
- Use cleaning and disinfectant agents according to the manufacturer's instructions.
- For patients with a diagnosis or suspicion of Creutzfeldt-Jakob disease (CJD) or variants, national regulations regarding disinfection and sterilization of products must be followed.
- Handle with validated methods and calibrated autoclaves and machines at all times.
- Prefer mechanical cleaning cleaning instead of manual cleaning as it gives more reliable results.
- Handle with care. Do not drop or apply excessive forces to the products during reprocessing.

6.1 Cleaning and Maintenance Products / Tools

warning The use of metal brushes and metal scourers causes corrosion!

warning The use of chlorine solutions (saline solutions, bleaching agent, ringer solution) causes pitting corrosion or stress corrosion cracking!

warning Use only process chemicals suitable for anodized aluminum and stainless steel.

- Long-term processes with disinfectants based on alkylamine active ingredient in silicone elastomers (silicone-coated products) may cause hardening.
- Use cleaning agents with approved efficacy (VAH or FDA approval or CE mark)
- Cleaners with a disinfectant effect (manual cleaning) must be in accordance with EN 14885 or equivalent national guidelines.
- Use soft, lint-free cloths, soft bristle brushes and compressed air guns.
- Do not use aldehyde-based cleaning products due to their protein-fixing properties.
- Cleaning solutions should be changed on daily basis. If the solution gets contaminated visibly, it can be changed earlier.

6.2 Water Qualifications

- **Drinking (Tap) Water***: Uses for the first wash and intermediate rinse.
 - Drinking water quality according to EC Directive 98/83/EC or AAMI TIR34
- **Softened water***: Water hardening substances (calcium and magnesium cations) are reduced.
- **Demineralised water***: Mineral substances largely removed by reverse osmosis or cation & anion exchangers or electrode ionization or distillation.

- For steam sterilization and the final rinsing of cleaning processing, limit values for feed Water quality is specified in EN 285 and ISO 17665.

⚠ Non-conforming water properties causes:

- Scaling and lime deposits,
- Pitting corrosion
- Brown-red deposits
- White-grey, colored surface changes
- Spotting and rusts.

7. CLEANING AND DISINFECTION INSTRUCTIONS

Preparation for Cleaning / Pre-cleaning

⚠ warning **Do not exceed** 25 °C for pre-cleaning temperature!

- Remove the lid from the container and separate all detachable components such as filter/filter retainer, basket, indicator and labels.
- Rinse all parts under cold tap water to remove coarse contaminations and visible soils.

Manual Cleaning / Disinfection

⚠ warning **Do not exceed** the maximum permitted cleaning temperature of 45 °C. Otherwise, protein denaturation may occur.

⚠ warning Use a non-protein fixing cleaning solution that **does not contain alcohol, aldehyde, phenol and QUAT** with a disinfection effect according to the manufacturer's instructions.

	Process Step	Water Quality	Heat	Time
1	Cleaning/disinfecting[1]	Drinking Water*	Room Temp.	2 min
2	First Rinsing	Drinking Water*	1 – 16° C	1 min
3	Disinfecting[1]	Drinking Water*	Room Temp.	10 min
4	Final Rinsing	Demineralized Water*	Room Temp.	2 min
5	Drying	-	Room Temp.	-

1. Pour a 25 ml dose the enzymatic enzymatic / neutral pH detergent solution in 5 liters of water (%0,5 dilution). Immerse the container completely in the solution and shake it gently to avoid air bubbles. Actuate all hinges and movable parts to ensure contact of the solution with all surfaces. Using a soft-bristled brush, remove all visible soil and debris while soaking the container.

2. Remove the containers from the solution and rinse under cold tap water for a minimum of 1 minute. Actuate all hinges and moveable parts while rinsing. Thoroughly and aggressively flush hard to clean areas.

3. Pour a 125 ml disinfectant solution dose in 5 liters of water (%2,5 dilution). Immerse fully in the disinfectant solution and allow container to soak for a minimum of 10 minutes. Actuate all hinges and movable parts to ensure contact of the solution with all surfaces. Using a soft-bristled brush, remove all visible soil and debris while soaking the container. Particular attention should be given to hinges, joints, and hard to clean areas.

4. Rinse with running demineralized water for at least 2 minutes.

5. Dry the products with a clean, lint-free cloth or compressed air.

[1] Anios - Aniosyme DD1 Solution was used to verify cleaning.

- Check for cleanliness visually at the end of the cleaning process especially the hinges, joints, etc. and hard to clean areas. Repeat the cleaning process if necessary.
- Apply maintenance procedures (see **Maintenance & Care**)
- Control visually and functioning of the containers according to the "**Inspection & Check**" part of this document.

Mechanical Cleaning and Disinfecting

⚠ warning Use low alkaline cleaning agents (8.0-10.5 pH). Do not use high alkaline value (>11.0 pH) cleaning agents.

⚠ warning Follow the washer-disinfector manufacturer's instructions. Use a washer-disinfector validated according to EN ISO 15883.

- Place the containers to the washer-difinfector in an open position and separated from their components.
- Arrange the position of the containers in the washer-disinfector to prevent mechanical damage through contact.
- Avoid overloading to make sure the surface of the containers effectively gets in contact with cleaning solution.
- With regard to the process chemicals used, consider the manufacturer's instructions concerning dosing, concentration, temperature, and exposure time.
- To prevent corrosion, unload the washer-disinfector immediately when the process steps done.

	Process Steps	Water Quality	Heat	Time
1	Pre-wash	Softened water*	<25 °C	2 min
2	Cleaning	Demineralized Water*	45 – 55 °C	5 min
3	First rinse / Neutralization	Demineralized Water*	>10 °C	2 min
4	Intermediate rinse	Demineralized Water*	>10 °C	2 min
5	Thermal disinfection / Final rinse	Demineralized Water*	90 °C	5 min
6	Drying	-	-	***

- Check for cleanliness visually at the end of the cleaning process especially the hinges, joints, etc. and hard to clean areas. Repeat the cleaning process if necessary.
- Apply maintenance procedures (see **Maintenance & Care**)
- Control visually and functioning of the containers according to the "**Inspection & Check**" part of this document.

***** See **"Water Qualifications"** section of this document for water properties.

****** If a high alkaline cleaning solution to be used, neutralization may be required.

******* Drying times vary depending on the loading capacity and must be observed according to the manufacturer's instructions.

Maintenance & Care

⚠ warning Ensure to apply care product (lubricant) to metal sliding surfaces, moving parts, joints prior to **"Inspection & Check"**.

⚠ warning Raw materials used in care products must comply with pharmacopeias such as USP/NF, Ph. Eur, and national and international cosmetics regulations. **Do not use** silicone oil-based products!

- Ensure that the products are cleaned effectively before maintenance. Visually and physically inspect the products for blood and protein soils.

Care products must be:

- Paraffin/white oil based,
- Biocompatible,
- Suitable for steam sterilization and vapor-permeable.

- Allow the products to cool down to room temperature.
- Apply the care product directly to joints, hinges and friction surfaces.
- Actuate moving mechanisms a few times to spread the care product evenly.
- Wipe the excessive care product with a lint-free cloth.

⚠ warning Failure to apply sufficient care products to moving parts causes friction-related damage and corrosion!

⚠ warning The silicone gaskets should not be treated with maintenance oil and solvents!

Inspection & Check

Verify the following:

- Check the laser marking of the product. It must be easily readable visually or scannable.
- The durability of the silicone gaskets is up to 500 sterilization cycles. Inspect the gasket to ensure that it is free of cracks and tears and that it is properly seated in its retaining groove. The silicone gasket must always be checked and replaced if required. Sealing should be checked before each use.
- There must not be any crack, fracture, wearing, burr, discoloration, spotting, rust or corrosion on the surface of the product.
- There must not be any loose, bent, broken, cracked, worn, or fractured parts.
- There must not be any discoloration or deterioration on the coating integrity of the surface, with anodized and colored products!
- Moving parts must be able to move freely without getting stuck or rubbing.
- After making sure that all parts perform their functions, assemble the components and check their compatibility with each other.

warning

If any of the above controls are not confirmed, **do not use** the product, contact with Bahadır Medical Devices (manufacturer). Damaged products **must only** be repaired by the manufacturer.

8. STERILIZATION INSTRUCTIONS

Sterilization Method

- The steam sterilizer must comply with the EN 285 standard and be verified according to the EN ISO 17665 standard.
- Ensure that the maximum permitted load according to the manufacturer's instructions is not exceeded.
- For steam sterilization, limit values for feed water quality as specified in EN 285 and ISO 17665.
- Ethylene oxide or other sterilization methods are not recommended.

Sterilization Method	Temperature	Holding Time	Drying Time
3x pre-vacuum steam sterilization	134 °C[1,3]	5 minutes[2,3]	20 minutes

[1] Validated exposure time required to achieve a 10-6 sterility assurance level (SAL)

[2] Validated exposure temperature required to achieve a 10-6 sterility assurance level (SAL).

[3] Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in this table.

Placing & Loading

warning

The total weight of a load of a container **must not exceed** the following loading sizes according to DIN EN 868-8. Otherwise, sterility assurance cannot be guaranteed!

warning

The containers **must not have** any additional packaging or be wrapped externally!

- According to EN 868-8 and DIN 58953-9, loading weight with instruments into a full size sterilization container (including basket/tray), should not exceed 10 kg due to ergonomic reasons and avoid excess condensation.
- After stacking the surgical instruments in the trays, place them in the container.
- In case of loading with textiles, pay attention that the pieces of laundry or folded textiles are positioned vertically and not be tight. DIN 58953-9
- Ensure that loaded materials do not prevent the working of the filter by keeping the required distance between the load and the lid.
- Perform the locking by placing the lid completely parallel to the container. Make sure the lid is fully seated on the container.
- Heavy and bigger containers should be positioned at the bottom of the sterilization chamber.
- When sterilizing, containers and textiles together, containers should be placed below absorbent, wrapped materials.
- The stacked containers on top of each other, must not exceed 45 cm in order to ensure effective air removal and steam penetration. Stacking is only recommended for fractionated vacuum method. The manufacturer's instructions must be followed.
- Always load containers by carrying them by their handles.
- After sterilization, allow the containers to cool down at room temperature.
- In accordance with national legislation, the loading limit may differ from the listed below.

8.1 Container Loading Capacities

Types	Dimensions (mm)	Recommended Instrument Loading (kg)	Recommended Textile Loading (kg)
Standard 1/1 Full Size	585x280x90	3,6	2,9
	585x280x120	4,8	3,8
	585x280x140	5,6	4,5
	585x280x190	7,6	6,1
	585x280x250	10,0	8,0
Standard 3/4 Middle Size	465x280x90	2,9	2,3
	465x280x120	3,9	3,1
	465x280x140	4,5	3,6
	465x280x190	6,1	4,9
	465x280x250	8,0	6,4
Standard 1/2 Half Size	295x280x90	1,8	1,4
	295x280x120	2,4	1,9
	295x280x140	2,8	2,2
	295x280x190	3,8	3,0
	295x280x250	5,0	4,0
Standard Wide Body	600x400x110	6,4	5,1
	600x400x170	9,9	7,9
Standard Long Body	710x280x190	9,2	7,4
Dental	310x190x80	1,3	1,0
	310x190x100	1,6	1,3
	310x190x130	2,1	1,7
Flat	300x280x45	1,0	0,8
	300x280x75	1,7	1,4
Mini	310x130x40	0,4	0,3
	310x130x70	0,7	0,6
	310x130x100	1,0	0,8
Mini Dental	210x155x50	0,4	0,3
	210x155x60	0,5	0,4

Storage

- Containers **must not be stored** together with chemicals!
- Sterile, packaged products should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.
- The containers maintain their sterility for 12 months under appropriate storage conditions. Tested in accordance with EN ISO 11607-1.
- Please refer to DIN 58953-8 (Logistics of sterile medical devices) for the storage period of medical devices in sterilization containers. The storage time usually depends on the storage conditions and must be determined by the responsible hygiene specialists.

Storage conditions:

- Temperature: 15-26°C
- Humidity: 30 – 50%
- Air pressure: 500-1060 hPa

Service Life

- Containers can provide service as long as they are regularly maintained, used according to their intended purpose, not damaged, and storage and maintenance conditions are complied with.
- 500 cycles are considered a minimum as a service life for containers and silicone seals can be used throughout the service life.
- Products that cannot pass functional and visual inspections should not be used.

warning

Products **must be sterilized** before being sent to the manufacturer or the manufacturer's technical service.

warning

In case of any serious incident related to the device, Bahadır Tıbbi Aletler AŞ (manufacturer) or local authorities should be reported.

Manufacturer Contact Information

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